



ES3M

**Manual Pure Tone
Audiometer**

Short-Form User's Guide

Rev 1.03




Micro Audiometrics
C O R P O R A T I O N

KEY USAGE

ES3M models have a yellow keypad/overlay with the following keys:

- ON** – Turns the instrument on
- MENU** – Menu access / Move up one level in Menu structure
- TALK** – Enter talk over mode
- Hz▲** – Increase frequency / Move up in menu list
- Hz▼** – Decrease frequency / Move down in menu list
- dB▲** – Increase level / Move up in menu list
- dB▼** – Decrease level / Move down in menu list
- L** – Present signal to Left ear / Select menu item
- R** – Present signal to Right ear / Select menu item
- L + R** – Present signal to Both ears

Press the  logo key to display the ES3M Information screen.
Press any key to exit the information display screen.

Key presses appear in { }; e.g. {**MENU**} means **press** the MENU key.

Use {**Hz▲**} / {**Hz▼**} or {**dB▲**} / {**dB▼**} for menu navigation.

Use {**L**} or {**R**} to select a menu item.

Menu sequences are represented as **←MENU1 ←MENU2 ←MENU3**.
The “←” means “scroll to Menu Item and press the **L** or **R** key”.

QUICK CHECKS

At power-up, the ES3 logo should be displayed and the backlight should be on.
After a short pause, the ES3M should display the **manual audiometry** screen and should respond to keypad control.

Keypad operation can be tested by pressing keys to see if the appropriate response occurs.
When in menu mode, the speaker should produce audible key press ‘ticks’.

Set frequency to 1000 Hz and level to maximum.
Press **L** or **R** and a tone (pulsed or continuous) should be heard at a comfortable loudness level at the correct ear.

Press **dB▲** and **dB▼** and present a signal to verify that signal loudness increases or decreases, respectively.

Press **Hz▲** and **Hz▼** and present a signal to verify that signal frequency increases or decreases, respectively.

MENU TREE (FROM TOP LEVEL MENU)

Audiometry Setup <ul style="list-style-type: none"> . Frequencies to Test . Starting Frequency . Levels to Test . Starting Level . Tone Mode . Lock or Unlock Settings . Reset Settings 	Audiometry options <ul style="list-style-type: none"> . Select from list . Select from list . Select from list . Select from list . Pulsed / Continuous . Password protected . Yes / No
Advanced Settings <ul style="list-style-type: none"> . A/C Power Settings <ul style="list-style-type: none"> . Backlight . Power Down . Bat. Power Settings <ul style="list-style-type: none"> . Backlight . Power Down . Adjust Contrast . Adjust Brightness . Key Volume . Perform Calibration <ul style="list-style-type: none"> . 70 dB HTL . 85 dB SPL . Calibration Date . Audiometer Test . Send Data 	Advanced setup options <ul style="list-style-type: none"> . When operating on A/C or USB power <ul style="list-style-type: none"> . Never / 30 sec / 1, 2, 5 min . Never / 1, 5, 15, 30 min / 1 hour . When operating on battery power <ul style="list-style-type: none"> . 5, 10, 20, 30 sec / 1 min . 15, 30 sec / 1, 2, 5 min . Adjust as desired . Adjust as desired . Low / Medium / High . Password protected <ul style="list-style-type: none"> . Set target level to 70 dB HTL . Set target level to 85 dB SPL . Enter Calibration Date . Tone “On” calibration test mode . Transmit cal data via serial port
Turn Off	Turn instrument off

Earscan® is a registered trademark of Micro Audiometrics Corporation.

Information in this manual is believed to be accurate; however, Micro Audiometrics Corporation assumes no responsibility for its use. Content is subject to change without notice. No part of this manual may be reproduced without written permission.