

**Earscan<sup>®</sup>**

**ES3**

Pure Tone  
Automatic Threshold Audiometer

**User's Manual**

Rev 1.13



**Micro Audiometrics**  
C O R P O R A T I O N

## **Precautionary Notes: Earscan® 3 Audiometer**

- 1. The Earscan® 3 is designed for use with (4) alkaline 1.5 volt AA size batteries.**
- 2. When installing batteries in the Earscan® 3, you must observe the correct polarity. If one or more cells are installed with reversed polarity, the instrument will not operate in the battery-powered mode, and the audiometer may be damaged.**
- 3. Use only the factory-supplied, medical grade power adapter (wall cube) with the Earscan® 3. Use of any other power adapter will invalidate the warranty and may result in damage to the audiometer.**
- 4. The Earscan® 3 is capable of producing signal levels greater than 100 dB SPL. Prolonged exposure to signals at these levels can result in temporary threshold shift or permanent hearing loss. When testing in manual mode at high HTL settings, limit signal duration to no more than 1 second.**
- 5. If the Earscan® 3 has only external power available (low batteries or no batteries installed), be aware that unplugging the mini DIN connector will result in power loss and potential data loss. If it is necessary to disconnect the mini DIN in this case (e.g., to switch from USB cable to printer cable), it is advisable to turn the instrument off, make the cable change, then turn the instrument back on.**
- 6. The Earscan® 3 does not implement handshaking protocol on the serial port. This means that the ES3 will not detect an unconnected serial port or a 'printer not ready' condition.**

# Table of Contents

<b>Introduction .....</b>	<b>5</b>
Manual Conventions .....	5
Keypad .....	6
Menu Navigation .....	6
Function Keys .....	7
<b>Setup .....</b>	<b>8</b>
Cable Connections .....	8
Power Up .....	8
Power Options .....	9
Beep Volumes .....	9
Talk Over .....	10
Time and Date .....	10
<b>Main Menu.....</b>	<b>11</b>
<b>Manual Audiometry .....</b>	<b>13</b>
Manual Audiometry Key Functions .....	14
Testing Procedure .....	15
Manual Audiometry Options .....	15
<i>Frequency Selection</i> .....	15
<i>Other Options</i> .....	16
<b>Automatic Threshold Audiometry .....</b>	<b>17</b>
Automatic Threshold Function Key .....	17
Test Procedure .....	18
Test Protocol .....	18
Automatic Threshold Test Error Conditions .....	19
Automatic Threshold Audiometry Options .....	19
<i>Frequency Selection</i> .....	19
<i>Frequency Presentation Order</i> .....	19
<i>Other Options</i> .....	20
<b>Test Results .....</b>	<b>22</b>
Displaying Results on the LCD .....	22
<i>Display Results Key Functions</i> .....	23
Quantifying Hearing Loss .....	23

**Serial Communications ..... 25**  
    Transmit Options ..... 25

**Screen and Clock Options ..... 26**  
    Screen Properties ..... 26  
        *Power-Up Logo* ..... 27  
    Real-Time Clock Settings ..... 27

**Lock Settings ..... 28**

**Reset Settings ..... 28**

**Calibration ..... 30**

**Technical Specifications ..... 31**

**Functional ‘Quick Checks’ ..... 32**

**Trouble Shooting Guide ..... 33**

**External Interface ..... 34**  
    Connectors ..... 34  
        *Multi-purpose Connector* ..... 34  
        *Headset Connector* ..... 35  
        *DPU-414 Printer Switch Settings* ..... 35

**Earscan 3 Information ..... 36**

**Warranty ..... 37**

**Contact Information ..... 38**

## List of Tables

Table 1. Power Settings .....	9
Table 2. Manual Audiometry Key Functions .....	14
Table 3. Other Manual Audiometry Options .....	16
Table 4. Automatic Threshold Audiometry Errors .....	19
Table 5. Basic Automatic Threshold Test Options .....	20
Table 6. Advanced Automatic Threshold Test Options .....	21
Table 7. Display Results Key Functions .....	23
Table 8. Scale of Hearing Loss .....	23
Table 9. Data Transmit Options .....	25

## List of Figures

Figure 1. ES3 Keypad .....	6
Figure 2. Mapping Function Keys to LCD Text .....	7
Figure 3. Manual Audiometry Display .....	13
Figure 4. Display of Threshold .....	14
Figure 5. Automatic Threshold Display .....	17
Figure 6. Audiometric Test Results Display .....	22
Figure 7. Information Display .....	36

## INTRODUCTION


The Earscan 3 model ES3 is a pure tone, full frequency, air conduction audiometer. ES3 audiometers are user customizable and can be battery, AC or USB powered. It provides manual/automatic threshold test capability with talk over. The ES3 is optimized for use in hearing conservation programs. It can be PC-controlled via optional software.

## Manual Conventions

The following conventions are used to indicate interactions with ES3:

KEY Presses appear in { }; e.g. {▲} means press the ▲ key.  
Basic menu navigation is done using {▲} / {▼} and {◀} / {▶}.

Use the {○} key to present tones, select menu items, or select list items.  
It's the navigation key in the center of the array:




Menu sequences appear as ◀Menu1 ▶Menu2 ▶Menu3. The “◀” symbol is used to indicate “scroll to Menu Item and press the ○ key, or press the number key corresponding to the Menu Item”.

# Keypad

ES3 audiometers have a blue keypad/overlay with 10 keys, labeled 1 through 0 (see Figure 1).

Figure 1. ES3 Keypad



## Menu Navigation

The ES3 menu system is designed for intuitive use. Navigating through menus is as simple as using the {▲}, {▼}, {▶}, or {◀} keys to highlight a selection and then pressing {O} to make a selection. The {Menu} key is used to exit test mode and enter the menu system, and also to move from 'lower' to 'higher' menu screens. Pressing {Menu} while in manual audiometry mode enters the 'top level' of the menu system, and if the user has navigated 'down' into the menu structure, each {Menu} press will return to the next 'higher' menu level until the 'top level' menu is reached. Menu sequences are shown assuming that the user begins at the 'top level' menu.


Note: The 'top level' menu can always be reached by repeatedly pressing {Menu}.

## Function Keys

At times, 'actions' or 'toggles' will be displayed on the bottom line of the LCD. The keys just below the display are then used as 'function' keys to perform an action (A), or to toggle between two possible parameter values (a/b).

**Figure 2** illustrates how keys are mapped to actions or toggles shown on the LCD.

**Figure 2. Mapping Function Keys to LCD Text**



In this illustrative example, there are four 'Actions' (Auto, Send, Clr, and Talk) assigned to 'Function' keys 1, 3, 4, and 5. Press **{1}** to perform action Auto, **{3}** to perform the action Send, **{4}** to perform action Clr, and **{5}** to perform action Talk. The 'L/R' shown above {2} indicates an L/R toggle with 'R' currently selected. Pressing {2} will toggle between the two possible states, and the state selected will be highlighted (e.g., **R**).



## SETUP

### Cable Connections

ES3 audiometers can be powered by 4 AA Alkaline batteries (battery compartment is accessible from the back of the instrument). Optionally, power from an AC adapter or USB bus power can be supplied via the 6-pin multi-purpose “mini DIN” connector on the top of the instrument. The AC adapter may be plugged directly into the ES3 mini DIN socket if serial communication is not required, or it may be plugged into the power adapter receptacle at the end of a printer cable or optional computer interface cable (RS-232). The headset and response button cable assembly plugs into the 15-pin connector on the top of the instrument.

Booth wall and headphone adapter cables are provided for installations that require patching through a booth. The booth adapter cable replaces the ES3 headphone/response button cable and terminates in three 1/4" phone plugs compatible with standard booth patch panels. The headphone adapter cable provides 1/4" mono phone plugs for TDH-39 headphones.

### Power Up

Press the **{On}** key, and the ES3 logo will scroll onto the screen unless scrolling has been disabled (see [Power Up Logo](#)), or the instrument was last powered down due to an inactivity time-out. When power-up initialization is completed, enter **{1}** to begin a new test or enter **{2}** to continue a previous test automatically saved by the ES3 at its last power down. If “New test” is pressed all current test results will be lost.

If you choose to begin a new test, the ES3 will prompt you to enter the Patient ID and an Operator ID. The manual audiometry screen will then be displayed and the instrument is ready to begin testing.

**Note:** When the ES3 powers down due to an inactivity timeout or in response to **Turn Off**, current test data is automatically saved to non-volatile memory. This prevents loss of data in case a test is interrupted. When the ES3S is powered back on, pending test data will be restored and the test can be continued, if desired, or a new test can be started.

# Power Options

Inactivity timeouts automatically turn off the LCD backlight and power after selectable periods of inactivity to conserve energy and extend battery life (see **Table 1**; ( ) indicates default setting).

**Note:** The backlight requires considerable battery power. Setting the backlight power-down interval longer than necessary will decrease battery life when the ES3 is battery powered (see Screen Properties).

**Table 1. Power Settings**

AC Backlight <b>▸Setup ▸Power ▸A/C ▸Backlight</b>	30 seconds 1 minute 2 minutes 5 minutes
AC Power Down <b>▸Setup ▸Power ▸A/C ▸Power Down</b>	Never 1 minute 5 minutes (15 minutes) 30 minutes 1 hour
Battery Backlight <b>▸Setup ▸Power ▸Battery ▸Backlight</b>	5 seconds 10 seconds (20 seconds) 30 seconds 1 minute
Battery Power Down <b>▸Setup ▸Power ▸Battery ▸Power Down</b>	15 seconds 30 seconds (1 minute) 2 minutes 5 minutes

# Beep Volumes

Key presses are silent when in test mode, but produce audible ‘beeps’ when in the menu system. The beep volume can be set via the menu sequence **▸Setup ▸Beep Volumes ▸Key Volume {Low/Medium/High}**. The ES3 also produces an alert sound when certain errors occur and to indicate the completion of an automatic test. The volume of the alert can be set via the menu sequence **▸Setup ▸Beep Volumes ▸Alert Volume {Low/Medium/High}**.

## Talk Over

Talk over mode may be entered by pressing {5} while in audiometric testing mode. There is no specific microphone 'opening'; sound arrives at the microphone through various openings in the cabinet (e.g., the slot in the side). Talk over volume may be adjusted with the {▲}, {▼}, {▶}, or {◀} keys. It is recommended that a normal speaking voice be used at a distance of 1 to 2 feet, and volume adjustment used to compensate for hearing status of the listener. Press {1}, {Menu}, or {O} to exit talk over.

## Time and Date

The ES3 has a battery-backed up real time clock to maintain time and date. The clock is set at the factory prior to shipment. See section [Real-Time Clock Settings](#) for details on how to check or set the time and date. The backup battery is a lithium rechargeable type and will be recharged whenever ES3 power is on.

**Note:** If the ES3 is not powered up for many months, the battery may become discharged. In this case, the time and date will need to be reset and the instrument should be left powered up for some time to recharge the clock battery. Normal instrument usage is sufficient to keep the battery fully charged.

# MAIN MENU

## 1-AUDIOMETRY

Enters manual test mode.

## 2-DISPLAY RESULTS

Displays test results on 1-2 screens.

## 3-DEMOGRAPHIC INFO

1-Patient ID

Up to 21 digits.

2-Operator ID

Up to 21 digits.

## 4-NEW TEST

Begins new test, overrides previous test.

**PRINT OR SEND PREVIOUS TEST BEFORE STARTING NEW TEST.**

## 5-SEND DATA

Transmits to printer or computer via serial port with supplied cables.

## 6-SETUP (Start Here)

Select options and save for all tests.

Can be changed when desired, resaved.

### 1-Audiometry Setup

### Audiometry options

1-Frequencies

Select from list

2-Min/Max Test Level

1-Minimum Test Level

-10 (default), -5, 0, 5, 10

2-Maximum Test Level

80, 85, 90, 95, 100 (default)

3-Manual Audiometry

Manual aud. options

1-Ear

**Left/Right**

2-Starting Level

10/15/**20**/25/30/35/40

3-Starting Frequency

Select from list (**1000**)

4-Tone Mode

**Pulsed**/Continuous

5-Increment Amount

**5**/10/20

6-Decrement Amount

5/**10**/20

7-Binaural Stimulus

**No**/Yes

4-Automatic Theshold

Automatic Test Options

1-Ear

**Left/Right**/Better Ear

2-Frequency Order

Select from list

3-Starting Level

10/15/20/25/**30**/35/40

4-Increment Amount

**5**/10/20/30

5-Decrement Amount

5/**10**/20/30

6-Gross Increment

5/10/**20**/30

7-Gross Decrement

5/10/**20**/30

8-Miscellaneous

Other automatic threshold options

1-Level Override

40 / 45 / 50 / 55 / 60 / 65 / 70

2-Faster Method

No / **Yes**

- 3-Stop On Inc.
- 4-Max Presentations
- 5-1kP Test Both Ears
- 6-Max False Responses
- 7-Ascend. Res. Required

- No / **Yes**
- 15 / **20** / 25 / 30 / 35 / 40
- No / **Yes**
- 1 / 2 / 3 / 4 / **5** / 10
- 1 / **2** / 3 / 4 / 5

## 2-Communications

- 1-Default Output
- 2-Baud Rate
- 3-Insert Linefeed
- 4-Auto Send Data

## 3-Date

- 1-Set Date Format
- 2-Set Date
- 3-Set Time Format
- 4-Set Time

## 4-Power

- 1-A/C
  - 1-Backlight
  - 2-Power Down
- 2-Battery
  - 1-Backlight
  - 2-Power Down

## 5-Display

- 1-Contrast
- 2-Brightness
- 3-Scroll Logo

## 6-Beep Volumes

- 1-Key Volume
- 2-Alert Volume

## 7-General

- 1-Calibration Dates
- 2-Perform Calibration
- 3-Headset for Testing
- 4-Earscan 3 Info.
- 5-Lock Settings

## 8-Reset Settings

## 7-TURN OFF

## Serial Communications Options

**Printer** / Computer

**9600** / 19.2 / 57.6 / 115.2

No / Yes

No / Yes

## Date Options

**MM/DD/YYYY** or DD/MM/YYYY

Enter date (use current date format)

24 Hour / 12 Hour

Enter time (use current time format)

## Power Handling Options

When operating on A/C or USB power

30 sec / 1, 2, 5 min

Never / 1, 5, **15**, 30 min / 1 hour

When operating on battery power

5, 10, **20**, 30 sec / 1 min

15, 30 sec / **1**, 2, 5 min

## Display Setting Options

Adjust as desired

Adjust as desired

**No** / Yes

## Volumes for key 'tick' and audible alerts

**Low** / Medium / High

Low / **Medium** / High

## Other Options and Information

Display last calibrated and due date

Password protected

Primary / Secondary

Displays information screen

Password protected

**Yes** / No

# MANUAL AUDIOMETRY

To insure that pending data is cleared, always begin a new test via the main menu sequence **↵New Test** to clear pending test data.

When the ES3 enters manual audiometry mode, the screen format shown in **Figure 3** will be displayed (assumes factory defaults for starting frequency and level).

**Figure 3. Manual Audiometry Display**

Frequency1000Hz

Level20dB

EarRight

Untested

AutoL / RSendClrTalk

1

2

3

4

5

6

7

8

9

0

Once a threshold has been obtained, the display will change to show the threshold value (the 'm' indicates threshold was obtained via manual testing; see **Figure 4**).

**Figure 4. Display of Threshold**

F r e q u e n c y	1 0 0 0	H z		
L e v e l	2 0 m	d B		
E a r	R i g h t			
A u t o	L / R	S e n d   C l r   T a l k		
1	2	3	4	5
6	7	8	9	0

The threshold may be cleared (Clr) by pressing key **{4}** (e.g., to retest a threshold). Only the currently-displayed threshold will be cleared.

The test ear, signal mode, starting frequency and starting level are all user selectable, so these values may differ from the ones shown. The example shown indicates that threshold of 20 dB has been obtained for the right ear at 1000 Hz.

## Manual Audiometry Key Functions

Table 2 lists the functions of ES3 keys that are active in the manual audiometry mode. The last column indicates the action of the function keys.

**Table 2. Manual Audiometry Key Functions**

{◀}/{▶}		Select next lower/higher enabled frequency (Hz)
{▲}/{▼}		Increase/decrease level by current step size (dB)
{1}	Auto	Exit manual test mode and enter automatic test mode
{2}	L/R	Toggle between Left and Right ear, or among Left, Right and Binaural if binaural testing is enabled
{3}	Send	Send the test data to computer or printer
{4}	Clr	Clear the current threshold
{5}	Talk	Exit test mode and enters talk over mode
{Menu}		Exits manual test mode and enters the menu system

## Testing Procedure

1. Instruct the patient to raise his/her hand or press the patient response button (if available) whenever a tone is heard.
2. Position the headphones over the patient's ears (Red over **Right** ear, Blue over **Left** ear).
3. Press **{2}** if necessary to toggle to the desired test ear.
4. Select the test frequency using the **{◀}** / **{▶}** keys.
5. Use the **{▲}** / **{▼}** keys to adjust the hearing level and the **{O}** key to present the stimulus to the selected ear. When testing with Continuous signals, stimulus duration should be about **1 second**; the minimum duration allowed is 200 ms.
6. Continue using the **{▲}** / **{▼}** and **{O}** to present stimuli consistent with the test paradigm being used (typically Hughson-Westlake) until threshold is determined.
7. Use the **{◀}** / **{▶}** keys to select the next frequency to test.
8. Repeat steps 4 through 7 until threshold has been determined for each frequency for the selected ear.
9. Change the test Ear by pressing **{2}**.
10. Repeat steps 4 through 8 until threshold has been determined for each frequency for the second ear.

## Manual Audiometry Options

### Frequency Selection

Note that the set of test frequencies selected for manual audiometry will also be used during automatic threshold testing. By default the frequency set includes 250, 500, 1000, 2000, 3000, 4000, 6000, and 8000 Hz, plus the 1k/P validity check. The 1k/P validity check is mainly used during an automatic threshold exam, but can also be tested during a manual exam if desired.

The frequency set to be used for manual or automatic threshold testing may be viewed or modified via the menu sequence **▸Setup ▸Audiometry Setup ▸Frequencies**. The display will show the list of available test frequencies and there will be a check mark (✓) next to the ones currently selected for testing. Use the **{▶}**, **{▲}**, **{▼}** and **{◀}** keys to move the highlight to a frequency to select or deselect for testing, and press **{O}** to toggle between 'selected' and 'deselected'. Repeat this process as necessary to select or deselect other frequencies.

Once selections are completed, press **{Menu}** to exit frequency selection mode and return to the menu system.



## Other Options

Additional manual audiometry options listed in **Table 3** may be accessed via the menu sequence **▢Setup ▢Audiometry Setup ▢Manual Audiometry**. To change an option, highlight the option to modify using the {▶}, {▲}, {▼} and {◀} keys, and press {O}. Use the {▶}, {▲}, {▼} and {◀} keys to select the desired setting for the option, and press {O}.

**Table 3. Other Manual Audiometry Options**

Ear	Ear to be tested first.	Left, (Right)
Starting Level	Initial level when new test is started	10, 15, (20), 25, 30, 35, 40
Starting Frequency	Initial frequency when entering manual test mode	125, 250, 500, (1000), 1500, 2000, 3000, 4000, 6000, 8000, 1k/P
Tone Mode	Stimulus mode	(Pulsed), Continuous
Increment Amount	Increment step size	(5), 10, 20
Decrement Amount	Decrement step size	5, (10), 20
Binaural Stimulus	Enable presentation to both ears simultaneously	(No), Yes

# AUTOMATIC THRESHOLD AUDIOMETRY

To insure that pending data is cleared, always begin a new test via the main menu sequence **↵New Test** to clear pending test data.

The automatic threshold test option performs a modified Hughson-Westlake test. **Figure 5** shows the display format used for Automatic Threshold tests. This example indicates that a threshold of 25 dB was registered in automatic test mode (indicated by “a”) for 3000 Hz at the left ear.

**Figure 5. Automatic Threshold Display**

F r e q u e n c y	3 0 0 0	H z		
L e v e l	2 5 a	d B		
E a r	L e f t			
M a n				
1	2	3	4	5
6	7	8	9	0

## Automatic Threshold Function Key

Pressing **{1}** (Man) during an automatic threshold exam will stop the test and return the ES3 to the manual audiometry mode. Note that exiting the test does not clear test data. The test may be re-entered and continued from the point at which it was interrupted.

## Test Procedure

1. Instruct the patient to press the response button whenever a tone is heard.
2. Position the headphones over the patient's ears.
3. From the manual audiometry screen, press key {1} to begin the automatic threshold exam.

## Test Protocol

1. The automatic test begins at the level and frequency indicated in the automatic threshold audiometry settings.

If the Faster Method setting is selected:

The stimulus increment size is set to the Gross Increment size in the automatic threshold audiometry settings. The increment size remains at this setting until the first response by the patient to a stimulus. After the first response, the increment step size is set to the Increment Amount in the automatic threshold settings. The stimulus decrement size is set to the Gross Decrement size in the automatic threshold audiometry settings. The decrement size will remain at this setting until the first negative response by the patient to a stimulus. After the first negative response, the decrement size is set to the Decrement Amount in the automatic threshold settings.

If the Faster Method setting is not selected:

The stimulus increment size is set to the Increment Amount and the decrement size is set to the Decrement Amount in the automatic threshold settings.

2. After each stimulus presentation, ES3 will either decrease the intensity level by the decrement size if the patient responded or increment the intensity level by the increment size if the patient did not respond.
3. The number of consecutive ascending positive responses at a given level is selectable in the automatic threshold audiometry settings. An ascending positive response is a positive response to a stimulus presentation that immediately follows a negative response to the previous stimulus presentation.
4. If the stimulus is presented at the maximum level for a given frequency and a negative response is received by the ES3, the threshold for the frequency is set to a No Response (NR) and the test continues with the next frequency selected for testing.
5. Testing for the next frequency in the test sequence will automatically start at an intensity level 10 dB HTL greater than the threshold for the previous frequency. However, the test level will not exceed the maximum starting level contained in the Level Override setting in the automatic threshold audiometry settings.

6. After testing at all the selected frequencies for the ear, the ES3 will examine the other ear to see if testing has been completed for each frequency. If the opposite ear has not been tested, ES3 will select the opposite ear, and begin testing in the manner described beginning with step 1. If both ears have been tested, the test is complete and ES3 will return to the manual audiometry mode.

## Automatic Threshold Test Error Conditions

If an error is encountered during an automatic threshold test, the ES3 will stop the test and indicate the error condition. The operator may re-instruct the patient via talk over mode and resume testing, or return to manual audiometry mode. **Table 4** lists the error conditions that the ES3 detects.

**Table 4. Automatic Threshold Audiometry Errors**

Button Held	Patient is holding down response button.
False Response	Patient pressed response button outside the valid response “window” or when no tone was presented.
Inconsistent Results	A valid threshold was not obtained. Testing stopped at the point where the error occurred.
Validity Error	The 1k and 1k/P thresholds differ by more than 5 dB.

## Automatic Threshold Audiometry Options

### Frequency Selection

Frequencies used during automatic threshold testing are the same as those for a manual test (see [Selecting Frequencies](#)).

### Frequency Presentation Order

Frequency presentation order may be changed via the menu sequence **▸Setup ▸Audiometry Setup ▸Auto Threshold ▸Frequency Order**. The default sequence is 1k/V, 500, 1000, 2000, 3000, 4000, 6000, and 8000 Hz. Only those frequencies selected for testing will be used.

To change presentation order, move the highlight with the {◀} / {▶} or {▲} / {▼} keys and press {O} to select a frequency. Move it to any position in the list using the {◀} / {▶} or {▲} / {▼} keys, and press the {O} key to place it at that position. When all frequencies are in the desired presentation order, press {Menu} to exit.

## Other Options

Basic automatic threshold test options (see **Table 5**) are selectable via the menu sequence **▸Setup ▸Audiometry Setup ▸Auto Threshold**.

**Table 5. Basic Automatic Threshold Test Options**

Ear	Starting ear for automatic threshold exam	Left, (Right), Better Ear
Frequency Order	Frequency presentation order	1k/P, 500, 1000, 2000, 3000, 4000, 6000, 8000
Starting Level	Starting level for automatic threshold test	10, 15, 20, 25, (30), 35, 40
Increment Amount	Level increment size	(5), 10, 20, 30
Decrement Amount	Level decrement size	5, (10), 20, 30
Gross Increment	Increment step size to use until positive response has been registered	5, 10, (20), 30
Gross Decrement	Decrement step size to use until negative response has been registered	5, 10, (20), 30

Advanced automatic threshold test options (see **Table 6**) are selectable via the menu sequence **▸Setup ▸Audiometry Setup ▸Auto Threshold ▸Miscellaneous**.

**Table 6. Advanced Automatic Threshold Test Options**

Level Override	Maximum starting level when changing frequency during automatic threshold exam.	40, 45, 50, 55, 60, (65), 70
Faster Method	Should gross increments and decrements be used?	No, (Yes)
Stop On Inc.	Should test be stopped if a reliable threshold cannot be determined?	No, (Yes)
Max Presentations	Maximum stimuli allowed per frequency before reporting 'inconsistent' error.	15, (20), 25, 30, 35, 40
1k/V Test Both Ears	Should 1k/Validity be tested for both ears?	No, (Yes)
Max False Responses	Maximum 'false responses' allowed before reporting an error.	1, 2, 3, 4, (5), 10
Ascending Responses Req	Number of consecutive ascending responses required to determine threshold.	1, (2), 3, 4, 5


# TEST RESULTS

Audiometric test results can be displayed on the LCD or transmitted to an attached PC or printer.

## Displaying Results on the LCD

Test results are displayed via the menu sequence **►Display Results**. Results obtained via automatic testing are displayed with an 'a', and results obtained via manual testing are displayed with an 'm'. **Figure 6** shows partial results obtained in manual test mode. The {◀}, {▶}, {▲}, or {▼} keys are used to scroll up or down through the results list.

**Figure 6. Audiometric Test Results Display**



## Display Results Key Functions

**Table 7** lists the valid keys when in the display test results mode.

**Table 7. Display Results Key Functions**

{◀}, {▶}	Scroll the display screen
{▲}, {▼}	
{Menu}	Exit display mode and return to the menu system

## Quantifying Hearing Loss

**Table 8** provides a general reference for converting threshold in decibels to degree of hearing loss.

**Table 8. Scale of Hearing Loss**

0 – 20 dB	Hearing within normal limits
25 – 40 dB	Slight to mild hearing loss
45 – 55 dB	Moderate hearing loss
60 – 70 dB	Moderately severe hearing loss
75 – 90 dB	Severe hearing loss
90 dB+	Profound hearing loss



(this page intentionally left blank)

## SERIAL COMMUNICATIONS

The ES3 can transmit data to a computer or to a printer. The default output device is selected using the menu sequence **▸Setup ▸Communications ▸Default Output**, and data is transmitted using the **▸Send Data** main menu option. Note that the ES3 implements **no handshaking protocol** for serial communications. This means that the ES3 cannot detect an unconnected serial port or a 'not ready' state. It is the user's responsibility to verify that a serial cable is connected and that the external device (printer or computer) is powered up and ready for communications before attempting to transmit data.

**NOTE:** While an Earscan 3 is attached to the PC using the supplied USB cable, it is automatically configured for output to a PC at a baud rate of 9600, and no setup changes are required.

### Transmit Options

**Table 9** lists optional settings associated with transmitting data.

**Table 9. Data Transmit Options**

Default Output	Output data destination	(Printer), Computer
Baud Rate	Serial port BAUD rate	(9600), 28.8k, 57.6k, 115.2k
Insert Linefeed	Insert linefeed after each line (for printer use)	No, (Yes)
Auto Send Data	Automatically send completed test to attached printer or PC	No, (Yes)

## SCREEN AND CLOCK OPTIONS

LCD brightness and contrast are user adjustable to allow optimizing the display for differing viewing conditions (e.g., ambient lighting or viewing angle). The contrast setting has no appreciable effect on battery life, but higher brightness settings require more power and will decrease battery life. The amount of time that the backlight remains on during periods of inactivity is user selectable, and should be set to as low a value as is convenient to preserve battery life.

The ES3 includes a real-time clock with rechargeable lithium battery backup to maintain time and date. Time can be set/displayed in either 12- or 24-hour format. Date can be set/displayed in dd/mm/yyyy or mm/dd/yyyy format. The battery is recharged when the ES3 is powered on. When fully charged it should maintain the real-time clock for at least a year, even if the instrument is not used. In the unlikely event that the clock battery should become discharged, the time and date would need to be reset and the instrument left on long enough to recharge the backup battery.

## Screen Properties

The LCD Screen is preset at the factory with typical contrast and brightness settings. Adjusting either setting may necessitate adjusting the other; e.g., increasing screen brightness may require changing contrast for optimum viewing.

LCD contrast can be adjusted via the menu sequence **▸Setup ▸Display ▸Contrast**. Use the {◀} / {▶} keys to adjust contrast to the desired setting. Press {O} to save the setting, or press {Menu} to exit without changing the contrast setting.

LCD brightness can be adjusted via the menu sequence **▸Setup ▸Display ▸Brightness**. Use the {◀} / {▶} keys to increase or decrease brightness. Press {O} to save the setting, or press {Menu} to exit without changing the brightness setting.

**Note:** A brighter backlight requires more battery power. Setting the backlight brightness to a value greater than needed will decrease battery life.

## Power-Up Logo

At power up, the ES3 logo can either scroll slowly onto the screen or just appear as a full-screen display. This option is set via the menu sequence **▸Setup ▸Display ▸Scroll Logo**. Select “Yes” for the scrolling display, or “No” for the full-screen display.

**Note:** If the ES3 powers down due to an inactivity timeout, the logo display will be skipped when the instrument is powered back up to minimize re-start time. Re-start will proceed directly to the manual audiometry screen.

## Real-Time Clock Settings

The date format is selected via the menu sequence **▸Setup ▸Date ▸Set Date Format**. The format may be set to **mm/dd/yyyy** or **dd/mm/yyyy**. The date may be set (or viewed) via the menu sequence **▸Setup ▸Date ▸Set Date**. A new date must be entered in the currently selected date format (i.e., dd/mm/yyyy or mm/dd/yyyy). Press {O} to save the date setting, or {Menu} to exit without changing the time setting.

The time format (12 or 24 hour) can be selected via the menu sequence **▸Setup ▸Date ▸Set Time Format**. The ES3 real time clock can be set (or time viewed) via the menu sequence **▸Setup ▸Date ▸Set Time**. A new time must be entered in the currently selected format (12- or 24-hour). Press {O} to save the time setting, or {Menu} to exit without changing the time setting.

## LOCK SETTINGS

All user selectable instrument settings may be “locked” if desired via the menu sequence **▸Setup ▸General ▸Lock Settings**. Enter a password and press **{O}** to lock all user settings at their current state.

**Caution!** Be sure to record or memorize the password used to lock instrument settings. It will be required to unlock settings.

To unlock settings, enter the menu sequence **▸Setup (password) ▸General ▸Unlock Settings** and enter the password that was used to lock the settings.

## RESET SETTINGS

All user selectable instrument settings may be returned to factory default values via the menu sequence **▸Setup ▸Reset Settings**. Calibration data will not be affected.

**Caution!** Applying ‘Reset Settings’ will return all user-selectable settings to factory defaults. User settings that differ from factory defaults will be lost.

(this page intentionally left blank)

## CALIBRATION

Micro Audiometrics recommends that audiometers be calibrated annually. About one month before the annual calibration due date, the ES3 will begin displaying a reminder at power up that the due date is approaching. This reminder will be cleared once the instrument is calibrated. If the calibration due date passes and the instrument has not been calibrated, a caution will be displayed regarding potential test data invalidation.

The last calibration date and calibration due dates may be viewed at any time via the menu sequence **▸Setup ▸General ▸Calibration Dates**. Please see the ES3 Calibration Guide for more detailed information regarding headset management and calibration.

# TECHNICAL SPECIFICATIONS

**ANSI S3.6 Type:** 4  
**Frequencies (Hz):** 125, 250, 500, 750, 1000, 1500  
**(+/- 1%)** 2000, 3000, 4000, 6000, 8000

## Levels

<b>TDH-39</b>	<b>125 Hz</b>	<b>-10 to 65</b>
<b>(60 Ω)</b>	<b>250 Hz</b>	<b>-10 to 90</b>
<b>Phones</b>	<b>500 thru 6000 Hz</b>	<b>-10 to 100</b>
<b>(+/- 1 dB)</b>	<b>8000 Hz</b>	<b>-10 to 90</b>

**Attenuation:** 5, 10, 20 db steps  
**Presentation:** Pulsed or Continuous  
**Repetition Rate:** Random intervals  
**Test Modes:** Manual, Automatic Screening / Threshold

**Display:** 128 x 64 Backlit LCD  
**Clock:** Real-time with lithium battery backup  
**Communications:** USB, RS-232  
**Headset:** TDH-39 60 Ω  
**Power:** 4 AA Alkaline Batteries, A/C adapter, USB Bus

**Standards Met:** ANSI S3.6-1996, ANSI/AAMI ES1:1993



## FUNCTIONAL ‘QUICK CHECKS’

The following checks can be used to verify instrument operation and to help narrow the focus for solving problems.

At power-up after a manual power-down, the ES3 logo should either scroll onto the LCD or appear quickly, depending on the scroll setting, and the backlight should be on. After an inactivity timeout power-down, the logo is bypassed. In either case, the ES3 should proceed to the “Select An Option:” screen and should respond to keypad control.

Keypad operation can be tested by pressing keys and verifying that the appropriate response occurs (e.g., pressing {▲} increases level, {▶} increases frequency, etc.). The speaker should produce audible ‘ticks’ when keys are pressed while in the menu system (keys are silent in test mode).

Next, check signal generation and keypad control. Set frequency to 1000 Hz, set level to 70 dB HTL, and press {O}. A tone (pulsed or continuous) should be heard at a comfortable loudness level. Verify that the signal is heard at the correct ear. Press {▲} {O} and {▼} {O} to verify that signal loudness increases or decreases, respectively. Press {2} {O} and verify that the signal is presented to the other ear.

Plug in a **response** button and verify that pushing the button causes the “ \* “ symbol to appear near the center of the display.

Enter the menu sequence **▸Setup ▸Date ▸Set Time** and verify that the time is correct and that the seconds count is incrementing.

Attach the **printer**, verify that printer output is selected via the menu sequence **▸Setup ▸Communications ▸Default Output**, and transmit data to the printer via the **▸Send Data** main menu option.

## TROUBLE SHOOTING GUIDE

Problem	Possible Solution
ES3 does not power up.	Check batteries or verify that wall cube is attached and plugged in.
Signal is missing or intermittent.	<ol style="list-style-type: none"> <li>1. Verify that headset/response button cable connector is securely attached and mounting screws are snug.</li> <li>2. Move or gently bend headset cable to see if problem “comes and goes” – if so, there may be a break in the wiring.</li> <li>3. Verify that screws holding the “fork” connectors at each earphone are snug.</li> </ol>
Subject response button does not work.	<ol style="list-style-type: none"> <li>1. Verify that headset/response button cable connector is securely attached and that mounting screws are snug.</li> <li>2. Disconnect response button plug from headset/response button cable, clean plug, and reconnect.</li> </ol>
Data is not properly sent to computer.	Reference: “Transmitting Data” <ol style="list-style-type: none"> <li>1. Verify that ‘output to computer’ is set.</li> <li>2. Verify that baud rate matches computer baud rate.</li> <li>3. Verify that interface cable is in place and computer is ready to receive data.</li> </ol>
Printer fails to print or prints incorrect characters.	Reference: “Transmitting Data” <ol style="list-style-type: none"> <li>1. Verify that ‘output to printer’ is set.</li> <li>2. Verify that printer is attached, power is “On”, and printer is ready to print.</li> </ol>
The Seiko printer does not properly feed paper.	Verify that paper is installed correctly and not jammed.
The Seiko printer appears to be printing, but nothing shows up on the paper.	Paper is in backwards. Thermal paper prints only on shiny side.

# EXTERNAL INTERFACE

## Connectors

The ES3 form factor necessitates a non-standard approach to headset and interface cable attachments. The headset attaches to a high-density 15-pin D-Subminiature connector, and is held in place with allen (hex) head screws to minimize the possibility of accidental mismatch of headset to instrument. All other cables are attached via a multi-purpose 6-pin mini DIN connector. The ES3 wall cube, printer/power cable, USB cable, or serial computer/power cable may be attached to the mini DIN connector. Connector pin outs are as follows:

### Multi-purpose Connector

6-pin mini DIN (pin configuration is shown as viewed from top of instrument)

Pin	Function	Pin	Function			
1	Ground	4	USB D+		6	5
2	+5 V In	5	RS-232 RX	4		3
3	RS-232 TX	6	USB D-		2	1

### Notes:

- 1) The wall cube connector uses pins 1 and 2 for power.
- 2) Serial cables use pins 3 and 5 for RS-232 and pins 1 and 2 for power (assuming wall cube is plugged into the cable's mini DIN 'pigtail').
- 3) The USB cable uses pins 1 and 2 for USB bus-power (requires about 120 ma) and pins 4 and 6 for USB data.

## Headset Connector

15-pin high density D-Sub (pin configuration is shown as viewed from top of instrument)

Pin	Function	Pin	Function
1	+5 V Out	9	Not Used
2	Lt Phn+	10	Not Used
3	Lt Phn-	11	Ground
4	Rt Phn+	12	Not Used
5	Rt Phn-	13	Not Used
6	RspButton	14	Not Used
7	Ground	15	Not Used
8	Not Used		

## DPU-414 Printer Switch Settings


The setup switches on the Seiko Instruments DPU-414 printer should be set as follows to ensure compatibility with Micro Audiometrics ES3 audiometers.

SW1 - 1 - OFF	SW2 - 1 - ON	SW3 - 1 - ON
2 - ON	2 - ON	2 - ON
3 - ON	3 - ON	3 - ON
4 - ON	4 - ON	4 - OFF
5 - ON	5 - ON	5 - OFF
6 - OFF	6 - ON	6 - ON
7 - ON	7 - ON	7 - ON
8 - ON	8 - OFF	8 - ON

## Earscan 3 Information

General Earscan 3 information is displayed via the menu sequence **▸Setup**  
**▸General ▸Earscan 3 Info.** (e.g., see Figure 7). Serial number and firmware revision are important when upgrades, updates, and/or service are being considered.

Figure 7. Information Display



# WARRANTY

## MICRO AUDIOMETRICS CORPORATION LIMITED WARRANTY

### EARSCAN 3 Pure Tone Audiometer

1. This is a "LIMITED WARRANTY" as defined in the Consumer Product Warranty and Federal Trade Commission Improvement Act. This WARRANTY gives you specific legal rights and you may also have other rights that vary from state to state.

2. Micro Audiometrics Corporation warrants this Earscan® 3 Pure Tone Audiometer to be free from defects in materials and workmanship for five (5) years and the headset and patient response button for one (1) year under normal use.

3. This WARRANTY does not cover items subject to normal wear and tear such as cables, earphone cushions, carrying cases, batteries, broken or marred cabinets, or any other accessories used in connection with this product, or consequential damages due to a defect in the product.

4. This WARRANTY does not apply to any product damaged by accident, misuse, tampering, alteration, abnormal condition of operation, carelessness, or if the products were altered or repair was attempted by anyone other than Micro Audiometrics Corporation or one of its Authorized Equipment Service Centers.

5. This WARRANTY applies only to the original customer, and only on units purchased and used solely within the United States and begins on the date of purchase. For your convenience keep the dated Invoice or Packing List as evidence of the purchase date.

Products not manufactured by Micro Audiometrics Corporation (noise reducing headphone enclosures, insert headphones, printers) are covered by their manufacturer's WARRANTY. Micro Audiometrics Corporation may, at its sole and exclusive option, repair or replace this product with either a new or like-new product provided that it has the functionality equal to the product replaced.

There are no obligations or liabilities on the part of Micro Audiometrics Corporation for consequential damages arising out of, or in connection with, the operation, use or performance of the product including, without limitation, with respect to loss of time, revenues or profits.

This WARRANTY does not cover transportation to and from the point of service. **"Loaner" equipment is not provided while service is being performed.**

\*\*\*\*\*

In the event of any claim of a defect covered by this WARRANTY, the customer should take the following steps:

- Contact Micro Audiometrics Corporation to discuss the nature of the claim.
- The audiometer and headset to be returned should be packed in the original shipping box. If not available choose an appropriate box with sufficient packing material to prevent damage during shipping.
- Return to Micro Audiometrics Corporation, 1901 Mason Ave, Suite 104, Daytona Beach, FL 32117 or to one of our Authorized Service Centers.

Micro Audiometrics Corporation disclaims all other warranties, expressed or implied, including any warranty of merchantability or for function of fitness for a particular purpose or application.

## **Contact Information**

**For additional information or assistance, contact your local distributor or contact Micro Audiometrics directly at:**

**Micro Audiometrics Corporation  
1901 Mason Ave, Suite 104  
Daytona Beach, FL 32117 USA**

**Phone: (386) 888-7878  
Toll-free: (866) 327-7226**

**Product and company information is available on the internet:  
[www.earscan.com](http://www.earscan.com)**

**For product information or inquiries, send email to:  
[sales@microaud.com](mailto:sales@microaud.com)**

**For product support or technical issues, send email to:  
[support@microaud.com](mailto:support@microaud.com)**



# **Micro Audiometrics**

C O R P O R A T I O N

1901 Mason Ave, Suite 104  
Daytona Beach, FL 32117  
Phone: (386) 888-7878  
**[www.microaud.com](http://www.microaud.com)**

**Earscan® is a registered trademark of Micro Audiometrics Corporation.**

**Information in this manual is believed to be accurate; however, Micro Audiometrics Corporation assumes no responsibility for its use. Content is subject to change without notice. No part of this manual may be reproduced without written permission.**

Copyright © 2004-2024 Micro Audiometrics Corporation.